

## Natus NatalCare LX — Advanced Incubator for Neonatal Intensive Care

The Natus NatalCare LX Infant Incubator provides an optimal thermal and developmental microenvironment for babies. Combining ease of use and accessibility for clinicians and families with economic value, the Natus NatalCare LX Infant Incubator meets the demands for different levels of acuity and the needs of various hospital departments.

The Natus NatalCare LX Infant Incubator can be used as a standard incubator or built-to-order with advanced options to meet the clinical and developmental needs of sick babies.

### RELIABLE THERMAL PERFORMANCE

- ▶ Double-wall acrylic canopy and doors minimize the radiant heat loss of the infant.
- ▶ Comfort Zone settings provide a range of recommended air temperatures based on the baby's birth age and weight.

### CONVENIENT PATIENT ACCESS & EASE OF USE

Ergonomic design provides the visibility and accessibility required for critical conditions, in addition to improving access to the baby during routine procedures.

- ▶ Large pull-out mattress tray, four access doors foldable 180°, 5 hand ports, and 8 IV ports offer quick and convenient access to the baby during procedures and therapies.
- ▶ Electronic tilt with automatic centering mechanism  $\pm 15^\circ$ , operating from the command module reduces the need to disturb the baby. Automatic centering feature returns the mattress to a flat position to help reduce errors during weighing, minimizing the need to re-weigh the baby.
- ▶ Easily positioned panel with large, easy to read LCD display provides optimal visibility of key parameters to the caregiver.

### EASY & EFFICIENT MAINTENANCE

The Natus NatalCare LX Infant Incubator is designed not only for ease of use but also for efficient cleaning. The incubator's smooth design and simple construction with only a few removable parts that are easy to disassemble and clean, minimizes equipment turnaround time.



# NatalCare LX

## infant incubator

The Natus NatalCare LX Infant Incubator provides an optimal thermal and developmental microenvironment for babies.

### All Natus NatalCare LX Incubators include the following:

- Mobile command module with LCD screen
- 2 skin temperature probes
- Trend graphs & comfort zone
- Electronic tilt with auto-centering
- Double-walled dome
- 4 access doors (2 double-walled) foldable 180°
- 5 oval ports

Part #	Item
KITNC-STLX0001	Natus NatalCare LX 120V
KITNC-STLX0003	Natus NatalCare LX with LIFT 120V — includes adjustable height lift option
KITNC-STLX0005	Natus NatalCare LX 120V NCA103 — includes servo-controlled humidity module
KITNC-STLX0007	Natus NatalCare LX LIFT 120V NCA103 — includes adjustable height lift option and servo-controlled humidity module

### Optional Modules

Part #	Item
NCA-102C	Oxygen servo-control module
NCA-103B/1	Humidity servo-control module 120V
NCA-105E	Electronic weighing scale module

### Accessories

Part #	Item
NCA-106C	Drawer
NCA-107	IV pole
NCA-108	Shelf
NCA-109	Right lateral door with iris port
NCA-109/1	Left lateral door with iris port
NCA-110	Antistatic wheels (4 with brake)
NCA-112	Short pole for accessories
NCA-113	Gel mattress
NCA-121	Type E oxygen cylinder holder
NCA-126	Flexible holder for ventilator tubing
NCA-128	Reusable water tank for NCA-103
NCA-131	Column for accessories
NCA-132/1	Universal auxiliary outlet 120V
SMC-NC01	Resuscitation system



### Supplies/Consumables

Part #	Item
DIR-06	Iris port sleeve (100/pk)
DIR-10	Skin temperature probe (reusable)
DIR-30	Disposable skin temperature probes (12/pk)
DIR-34	Skin probe protectors (120/pk)
DIR-36	Hand port sleeves (100/pk)
R18130	Servo humidity hose assembly
R18130/10	Servo humidity hose assemblies (10/pk)
TR-69	Air filter (4/pk)