

## OXYGEN SENSOR REPLACEMENT GUIDE



**C-1R**  
B39022-C1R



**C-2R**  
B39023-C2R



**R-13**  
B39636-R13



**R-15**  
C41131-R15



**R-17MED**  
C43690-R17MED



**R-22MED**  
C44611-R22MED



**R-23MED**  
C61339-R23MED



**R-24MED**  
C61340-R24MED



**R-26**  
C64335-R26MED



**R-29MED**  
C70592-R29MED



**R-30MED**  
B72938-R30MED



**R-31MED**  
B73460-R31MED



**R-32MED**  
B73463-R32MED



**T-1**  
A37016-T1



**T-2**  
A38070-T2



**T-4**  
A37017-T4



**T-7**  
A51327-T7

MANUFACTURER	MODEL/REF. NO.	SENSOR CLASS
Alphamed		R-17
Ametek	S-3A / 1 / ii	Not available
Analytical Industries	2000	Not available
Atom	0X-60	R-17
BCI		R-17
BMD		R-17
Bertucci		T-1 or T-2
BIO MS	Airox	R-29
BioMed	M-2 / M10, 48351	T-1 or T-2
Bio-Tek	X74223	T-1 or T-2
Biotest	S101	R-15
	C102-1	R-17
	A103	R-22
	406931	Not available
Bird	Sentry	R-17
BMD	10006	R-17
Burke & Burke	OOM104	R-26
	OOM103A	R-22
	OOM201	R-23
Caradyne	Criterion 40, 60	R-26
Catalyst Res. (See MSA)		
Ceramtec (See Maxtec)		
Cheiron		R-24
CIG Healthcare		T-4
CR		R-15
Criticare Systems		Not available
Critikon		Not available
Dameca	42035-46	R-17
Datascope	0600-00-0002	R-22
Datex Ohmeda	OT-102	Not available
	5250	R-15
	5250RGM	R-15
	4700 Oxicap	R-15
	ExcelMod.	R-15
	5100	R-15
	5120	R-15
	5125	R-15
	5400	R-15
	7800, 0237-2034-700	R-15
	7900, 6050-0004-110	R-30
	Modulus Excel	R-32
	Extended Life Sensor	R-32
Diversified Diag.	Xenematic	T-1 or T-2
Dp Medical	MVA 210	R-23
Drager	Babilog 8000	R-23
	S sensor	Not available
	Anemone	R-23
	Cicero	R-23
	Cato	Not available
	Evita	R-23
	Babilog Sa2	R-23
	Oxydig	R-23
	Oxycom	R-23
	PM8030	R-23
	I-8000	Not available
	8000IC	R-23
	Narkomed Series	R-31
Drager (N. America)	6803290 Catco	R-15
Emerson	3MV-OA	T-1
Engstrom (Datex)	Erica, Evita	R-17
Envitec	OxiQuantM	R-17
	OxiQuantRM	R-15
	OxiQuantS	R-22
F. Stephan		R-17
Fresenius	X74223	T-1
Hamilton	HM-01	T-1 or T-2
	HM-03	T-7
	HM-10	R-15
	HM-11	R-17
	HM-12	R-22
	HM-13	R-23
	HM-14	R-23
Henleys	5569	T-1
	5567	R-17
	5556	R-22
Hewlett Packard	HP15201A	R-15 or R-23
Heyer	Pocket O2	T-1 or T-2
Hill-Rom AS		R-22

MANUFACTURER	MODEL/REF. NO.	SENSOR CLASS
Hudson RCI / Ventronics	5500	R-13
	5525	C-1R or C-2R
	5528	C-1R or C-2R
	5529	C-1R or C-2R
	5550	C-1R or C-2R
	5552	C-1R or C-2R
	5556	R-15
	5557	R-23
	5558	T-7
	5565	C-1R or C-2R
	5566	R-22
	5567	R-17
	5568	T-1
	5569	C-1R or C-2R
	5570	T-4
	5575	T-4
	5577	R-13
	5584	R-13
	5590	R-13
	5595	R-13
	6477	R-13
Imed	810, 825	T-1 or T-2
Infrasonics		R-22
Inmed		T-1 or T-2
Ivac	1200, 1225	T-1 or T-2
Libra	OM500	T-1 or T-2
	OM600	T-4
Marquette / Hellige		R-22
Maxtec	OM-25AE & OM-25ME	R-29MED
	MaxO2 OM-25A	R-29IMED
Medigas		T-1 or T-2
Megamed	M1	C-1R or C-2R
	M2	T-1 or T-2
MSA	Mini-Ox I, II, III	R-17
	MinIA	Not available
Newport Medical	OM100	T-7
Nivaco		R-15
Norco Med		R-22
Novametrics		T-4
Ohmeda (See Datex)		
Oxyquip		T-1 or T-2
Oxitron	3710	R-23
Pacifictech	PT-10	R-15
	PT-14	R-23
Patient Tech		T-1 or T-2
PPG	250-2000-002 (SaraCalp)	R-22
PPG Hellige	5055055-001 (Irisa)	R-22
Prolab	OOM105	T-7
Puritan Bennet	7200	R-17
	7820	R-17
	4-020933-00	R-17
	PB850	R-17
Schoch		T-1 or T-2
Sechrist	650	T-7
Shock		T-1 or T-2
Siemens	Servo 300	R-24
	V710 300	R-24
	Servo 700	R-24
	Servo 900C	R-24
	Servo 900D	R-24
Spacelabs	90514	R-17
Sun Medical		T-1
Taema (CFPO)	Cesar YR010100	R-17
	Oxy2100	T-1
	Cesar (new)	R-22
Teledyne	TED-60T	T-7
	TED-160	T-1 or T-2
	TED-190	T-1 or T-2
	TED-191	T-7
	TED-200T	T-7
	TED-200T7	T-7
Toptronics	XT74223	T-1
Ventronics (See Hudson)		
Vickers	251	T-1 or T-2
VTI	100352	R-15
	100580	R-17
	10372	R-23
	VT1272	R-15
Waidray	MRI	T-7

NOTE: All sensors come with a one-year warranty against defects in material and workmanship.

### OXYGEN ANALYZERS AND MONITORS

TED 60-T      TED 191      TED 200-T7      Coming Soon

### TELEDYNE ANALYTICAL INSTRUMENTS

16830 Chestnut Street  
City of Industry, CA 91748, USA

TEL: 626-934-1500 FAX: 626-934-1651  
TOLL FREE: 888-789-8168

[www.teledyne-ai.com](http://www.teledyne-ai.com)